

Adjustable Tissue Tray

Part #: ATT-10

Price: \$175.00 USD

Size: 3.25" x 12" x 7.5"

Weight: 1 lb

The Adjustable Tissue Tray comes with one Subcuticular Suturing Model (SCS-10) and is compatible with BOSS. Its ergonomic design promotes proper hand placement for suturing practice.



Skill Development:

- » Basic and advanced suturing
- » Suturing open wounds under tension

Features and Benefits:

- » Tension device allows user to adjust pressure on suturing site
- » Tray has been designed to nest in the large open site on BOSS covering the alligator clips
- » Ergonomically designed for hand placement
- » Compatible with all of Simulab's tissue pads

Product Components:

- » Adjustable Tissue Tray
- » Subcuticular Suturing Model