

# Basic Open Surgical Skills (BOSS)

**Part #:** BOSS-10

**Price:** \$225.00 USD

**Size:** 17" h; 13" w; 4" d

**Weight:** 5 lbs

Basic Open Surgical Skills (BOSS) is a flexible platform for basic surgical skill development. It provides an easy to use platform for synthetic or cadaveric task modules and its ergonomic design creates a realistic tool for improving skills.

## Skill Development:

- » Basic and advanced suturing
- » Knot tying under tension
- » Knot tying at depth

## Features and Benefits:

- » Durable multi-use platform can be used with any combination of materials
- » Ergonomic design ensures correct hand positioning
- » Small, portable platform perfect for simulation lab or home use

## Product Components:

- » Tubing for knot tying under tension
- » Hook for knot tying
- » Cylindrical cup for knot tying and suturing at depth
- » Tissue tensioning device
- » Open site with flexible fixation system, including removable alligator clips and large 1" clips

## Compatible Products:

- » Student Tissue Pack (BTP-10)
- » Tissue Model (TSM-10)
- » Complex Tissue Model (TSC-10)
- » Large Intestine (LGI-10) or Small Intestine (SMI-10)
- » Tissue Suture Pad (TSP-10)

