

# FemoraLineMan<sup>®</sup> System

**Part #:** FLM-30

**Price:** \$1,230.00 USD

**Size:** 18" h; 27" w; 8.5"d (carrying case)

**Weight:** 27 lbs

FemoraLineMan offers an effective training solution for vascular access using the femoral vein. This trainer uses the same patented technology as the highly acclaimed TraumaMan System and allows medical professionals to train using real-time ultrasound guidance during catheter placement.

## Skill Development:

- » Vascular access using the femoral vein
- » Using ultrasound guidance during catheter placement
- » Identifying landmarks significant to the procedure

## Features and Benefits:

- » Anatomically correct body form with identifiable landmarks
- » Differentiates arterial and venous blood
- » Self sealing veins and skin for multiple cannulations
- » Replaceable tissue sets allow for multiple uses
- » Arterial pulse and proper landmarks are present to help avoid and detect errors
- » Ultrasound compatible for guidance during needle placement

## Product Components:

- |                                |                        |
|--------------------------------|------------------------|
| » FemoraLineMan Body Form      | » Arterial Fluid (red) |
| » Replaceable Tissue (FLMT-30) | » Vessel Filling Port  |
| » Pulse - Hand Pump            | » Fill Line            |
| » Venous Pressure Regulator    | » Carrying Case        |
| » Venous Fluid (blue)          | » User's Guide         |

