

FemoraLineMan Replacement Tissue


Part #: FLMT-30

Price: \$325.00 USD

Size: 8" h; 7" w; 3.25" d

Weight: 3 lbs

The replaceable tissue for FemoraLineMan (FLM-30) provides a realistic look and feel and comes prefilled with two colors of fluid to simulate a venous or arterial location. Each tissue lasts for approximately 20-25 full catheter insertions or 100 to 200 needle sticks and is ultrasound compatible for needle guidance.



Skill Development:

- » Central venous access using the femoral vein route
- » Using ultrasound guidance during catheter placement
- » Identifying landmarks significant to the procedure

Features and Benefits:

- » Differentiates arterial and venous blood to show positive or negative results
- » Self sealing veins and skin for multiple practices
- » Arterial pulse and proper landmarks are present to help avoid and detect errors
- » Ultrasound compatible for guidance during needle placement