

Large Intestine

Part #: LGI-10

Price: \$40.00 USD

Size: 12.5" (length)

Weight: .1 lbs

Realistic, reinforced single layered bowel section for suturing or anastomosis.

- » Outer diameter: 29mm
- » Inner diameter: 27mm
- » Wall thickness: 2mm

Skill Development:

- » Suturing
- » One and two layered Anastomosis
- » Cutting
- » Tissue Manipulation

Features and Benefits:

- » Realistic tissue simulation
- » Contains submucosa

Compatible Products:

- » LapTrainer (LTS-40)
- » BOSS (BOSS-10)
- » Laparotomy Model (LPM-10)

