

Nissen Fundoplication Model

Part #: LNF-10

Price: \$725.00 USD

Size: 10" h; 9" w; 4" d

Weight: 3 lbs

This procedural training module is designed to practice either Laparoscopic or Open Nissen Fundoplication. The model is designed to fit inside any of Simulab's abdominal trainers for laparoscopic or open procedural training (see list of compatible housing below). The anatomical, durable plastic bed includes visual landmarks and the soft tissue replaceable parts which allow each student to complete the procedure. The base of this product is also compatible with a replaceable tissue that includes a stomach with an esophagus and a replaceable lower esophageal sphincter (LNF-20).



Skill Development:

- » Laparoscopic Nissen Fundoplication
- » Open Nissen Fundoplication

Features and Benefits:

- » Model base with replaceable tissue is designed for multiple uses
- » Soft tissue base with anatomical landmarks
- » Procedures are completed on replaceable tissues

Compatible Laparoscopic Trainers:

- » LapTrainer with SimuVision (LTS-40)
- » Pop-up Trainer (SPT-20)
- » Pop-up Trainer with Scope Holder (SPT-30)
- » Torso Trainer (ST-10)