

# Laparotomy Model


**Part #:** LPM-10

**Price:** \$385.00 USD

**Size:** 15" h; 12" w; 6" d

**Weight:** 7 lbs

The laparotomy model simulates a partial abdomen. It includes a replaceable abdominal tissue cover that includes skin, subcutaneous fat, fascia, preperitoneal fat, and peritoneum. The base is water tight and can be filled with fluid. It also includes two alligator clips for securing internal tissue structures like a simulated intestine (LGI-10). Additional replaceable tissues (LPMT-10) are sold separately.



## **Skill Development:**

- » Midline abdominal incisions and closure
- » Bowel anastomosis
- » Surgery in the abdominal region

## **Features and Benefits:**

- » Lifelike synthetic tissue
- » Model base designed for multiple uses with replaceable tissue
- » Flexible model base can be used for a variety of procedures

## **Compatible Products:**

- » Large Intestine (LGI-10)
- » Replaceable Tissue (LPMT-10)