

LapTrainer with SimuVision


Part #: LTS-40

Price: \$1,795.00 USD

Size: 22" h; 18" w; 8" d (carrying case)

Weight: 17 lbs

This highly portable, innovative laparoscopic video trainer features SimuVision™ technology, which is a simulated laparoscope that leverages PC technology to create unsurpassed graphics. SimuVision technology integrated into our laparoscopic trainer creates an affordable and portable platform for product and surgical demonstrations or for teaching techniques and skill building tasks.



Skill Development:

- » Basic laparoscopic skills and dexterity
- » Improve eye-hand coordination
- » Practice transferring, dissection and suturing

Features and Benefits:

- » Boom mounted camera head can be adjusted to create 0°, 30°, or 45° scope angles
- » Excellent alternative to costly virtual trainers
- » Can be used with a wide variety of skill development exercises (see compatible products below)
- » Portable and affordable
- » Similar look and feel of an actual laparoscope
- » Sets up in less than five minutes and can be used with any PC or laptop
- » FLS compatible
- » Plugs into a USB 2.0 port - no extended power needed

Compatible Products:

- » LapTrainer Skills Set #1 (LTSSI-10)
- » Cholecystectomy Model (LC-10)
- » Nissen Fundoplication Model (LNF-10)

LapTrainer with SimuVision

Product Components:

- » LapTrainer Top (LTP-1001)
- » Set of Two Side Supports (LTP-1002)
- » White Vinyl Top (OMR-30)
- » Clear Vinyl Top (CMR-30)
- » Rear Support Bracket (LTP-1004)
- » LapTrainer Base (LTP-1005)
- » Laptop Platform (LTP-1006)
- » Laptop Support Bracket (LTP-1007)
- » Carrying Case (LTP-1008)
- » SimuVision Camera with integrated lighting (LTP-1009)
- » Camera Protector (LTP-1010)
- » SimuVision Set-up Disc (LTP-1015)
- » User's Guide (LTP-1016)

* Instruments and PC not included