



# LapTrainer with SimuVision Package

**Part #:** LTS-PACKAGE  
**Price:** \$1,925.00 USD  
**Size:** 22" h; 18" w; 8" d  
 (carrying case)  
**Weight:** 17 lbs

This highly portable, innovative laparoscopic video trainer features SimuVision™ technology, which is a simulated laparoscope that leverages PC technology to create unsurpassed graphics. SimuVision technology integrated into our laparoscopic trainer creates an affordable and portable platform for product and surgical demonstrations or for teaching techniques and skill building tasks.



This skill set includes four standardized exercises ranging from basic to more advanced in laparoscopic skill development. Perfect for use with the LapTrainer with SimuVision (LTS-40).

## Skill Development:

- » Basic laparoscopic skills and dexterity
- » Improve eye-hand coordination
- » Practice transferring, dissection and suturing

## Features and Benefits:

- » Boom mounted camera head can be adjusted to create 0°, 30°, or 45° scope angles
- » Excellent alternative to costly virtual trainers
- » Can be used with a wide variety of skill development exercises (see compatible products below)
- » Portable and affordable
- » Similar look and feel of an actual laparoscope
- » Sets up in less than five minutes and can be used with any PC or laptop
- » FLS compatible
- » Plugs into a USB 2.0 port - no extended power needed

# LapTrainer with SimuVision Package

## LapTrainer Product Components:

- » LapTrainer Top (LTP-1001)
- » Set of Two Side Supports (LTP-1002)
- » White Vinyl Top (OMR-30)
- » Clear Vinyl Top (CMR-30)
- » Rear Support Bracket (LTP-1004)
- » LapTrainer Base (LTP-1005)
- » Laptop Platform (LTP-1006)
- » Laptop Support Bracket (LTP-1007)
- » Carrying Case (LTP-1008)
- » SimuVision Camera with integrated lighting (LTP-1009)
- » Camera Protector (LTP-1010)
- » SimuVision Set-up Disc (LTP-1015)
- » User's Guide (LTP-1016)

## Skills Set Product Components:

- » LapTrainer Top (LTP-1001)
- » Set of Two Side Supports (LTP-1002)
- » White Vinyl Top (OMR-30)
- » Clear Vinyl Top (CMR-30)
- » Rear Support Bracket (LTP-1004)
- » LapTrainer Base (LTP-1005)
- » Laptop Platform (LTP-1006)
- » Laptop Support Bracket (LTP-1007)
- » Carrying Case (LTP-1008)
- » SimuVision Camera with integrated lighting (LTP-1009)
- » Camera Protector (LTP-1010)
- » SimuVision Set-up Disc (LTP-1015)
- » User's Guide (LTP-1016)

\* Instruments and PC not included