

LapTrainer Skills Set #1

Part #: LTSSI-10

Price: \$225.00 USD

Size: 13" h; 13" w; 4.5" d (carrying case)

Weight: 3 lbs

This skills set includes four standardized exercises ranging from basic to more advanced in laparoscopic skill development. Perfect for use with the LapTrainer with SimuVision (LTS-40).

Skill Development:

- » Laparoscopic coordination, technique, and precision

Features and Benefits:

- » The Pegboard develops eye-hand coordination and ambidexterity by requiring the user to lift each disk, transfer between hands and place disks on a designated peg
- » The Key Trainer develops 3D to 2D orientation skills, ambidexterity, and working through unnatural wrist rotations
- » The Tissue Suture Pad develops skills required for needle transferring, placement of a suture, and knot tying
- » The Multi-Purpose Pad which comes with a simulated intestine, allows the user to secure any soft tissue model in place using four corner mounted clips, and practice advanced suturing and immobilizing tissue skills

Product Components:

- » One pegboard
- » One key trainer (SKT-10)
- » Two tissue suture pads (TSP-10)
- » One multi-purpose pad with bowel section (LGI-10)
- » Carrying case

