

Knot Tying Board

Part #: SKB-10

Price: \$65.00 USD

Size: 10" h; 8" w; 2" d

Weight: 1 lb

The Knot Tying Board has a multipurpose platform for skill development. It can be used on its own or with Simulab's laparoscopic trainers. The board includes parallel knotting tubes that simulate a realistic tissue response. It also includes two suture pads that simulate different skin thicknesses.

Skill Development:

- » Knot tying
- » Knot tying under tension
- » Suturing

Features and Benefits:

- » Multiple tasks can be performed
- » Portable and lightweight

Compatible Products:

- » LapTrainer (LTS-30)
- » Pop-Up Trainer (SPT-20)
- » Pop-Up Trainer with Scope Holder (SPT-30)

