

Key Trainer


Part #: SKT-10

Price: \$25.00 USD

Size: 3" h; 3" w; 2" d

Weight: 1 lb

The key trainer provides a platform to practice hand-eye coordination. It can be used with any Simulab laparoscopic training models. It allows users to build skills and confidence by passing a specially designed key through a narrow slot. The base is easily Velcroed to any surface to allow trainees to vary the angle of the tasks and learn difficult tool angles with both hands.



Skill Development:

- » Hand-eye coordination
- » Dexterity with laparoscopic instruments

Features and Benefits:

- » Small and portable
- » Durable

Compatible Products:

- » LapTrainer (LTS-40)
- » Tower Trainer (TPT-10)
- » Torso Trainer (ST-10)