

# SonoMom<sup>®</sup> System

**Part #:** SMM-10

**Price:** \$6,950.00 USD\*

**Size:** 30" x 23" x 11"

**Weight:** 40 lbs

This ultrasound trainer provides an affordable platform for teaching students how to read obstetric ultrasound imaging. SonoMom's soft tissue body form includes internal and external landmarks and 144 unique probe locations providing a normal image for each window. The system also includes two simulated probes - transabdominal and transvaginal.



## Skill Development:

- » Read diagnostic ultrasound imaging
- » Identify external and internal landmarks for a diagnostic exam

## Features and Benefits:

- » Simulated ultrasound probes that plug into a PC
- » Software that allows any PC to simulate an ultrasound machine (PC/Laptop not included)
- » Findings include patients who are not pregnant, in varying stages of pregnancy, and patients with abnormal pregnancies.
- » Body form with embedded targets for each of the viewing windows, depending on the module

## Product Components:

- » SonoMom Body Form
- » Simulated Ultrasound Probe - Transabdominal
- » Probe Protector - Transabdominal
- » Simulated Ultrasound Probe - Transvaginal
- » Probe Protector - Transvaginal
- » SonoMom User's Guide
- » Carrying Case

## SonoMom Training Modules:

- » SonoMom First Trimester Module \$1,950

\*SonoMom must be purchased with one module