

SonoMan[®] System

Part #: SNM-10

Price: \$5,450.00 USD*

Size: 33" h; 23" w; 11" d

Weight: 33 lbs

The SonoMan System (SNM-10) includes a soft tissue body form with internal and external landmarks and a simulated probe. The torso has 258 unique probe locations providing a normal image for each window. The system's software platform allows users to expand SonoMan's capabilities with additional training modules.

Skill Development:

- » Read diagnostic ultrasound imaging
- » Identify external and internal landmarks for a variety of diagnostic exams

Features and Benefits:

- » Software that allows any PC to simulate an ultrasound machine (PC/Laptop not included)
- » Transverse and longitudinal viewing options
- » Video of actual patients
- » Normal to abnormal patient images

SonoMan System Product Components:

- » SonoMan Body Form
- » Simulab Ultrasound Probe
- » Carrying Case

SonoMan Training Modules:

» SonoMan AAA Training Module (SNMA-10)	\$1,250
» SonoMan ECHO Training Module (SNME-10)	\$1,950
» SonoMan eFAST Training Module (SNMEF-10)	\$1,050
» SonoMan FAST Training Module (SNMF-10)	\$1,950
» SonoMan Gallbladder Training Module (SNMG-10)	\$1,450
» SonoMan Renal Training Module (SNMR-10)	\$1,950
» FAST Module with SONORIO (TMF-10)	\$1,950

*SonoMan must be purchased with one module

