

Pop-Up Trainer


Part #: SPT-20

Price: \$195.00 USD

Size: 16" h; 16" w; 8.5" d

Weight: 6 lbs

The Pop-Up Trainer is a portable, lightweight box trainer with a soft cover and can be positioned in various heights and angles. This unit can be used with any combination of the five spacers provided to obtain the desired height or angle. The soft cover allows for numerous trocar placement sites.



Skill Development:

- » Laparoscopic training
- » Gynecological training

Features and Benefits:

- » Collapsible for easy storage
- » Affordable
- » Adjustable height
- » Personalized trocar placement
- » Great for reps or if storage space is limited

Compatible Products:

- » Cholecystectomy Model (LC-10)
- » Nissen Fundoplication Model (LNF-10)
- » Choledochotomy/CBDE Model (LCD-10)
- » Cholecystenterostomy Model (LCE-10)
- » Tissue Model (TSM-10)
- » Complex Tissue Model (TSC-10)
- » Tissue Suture Pad (TSP-10)