

Folding Suture Trainer

Part #: SST-20

Price: \$295.00 USD

Size: 13" h; 13" w; 4.5" d

Weight: 4 lbs

This small, portable trainer is perfect for laparoscopic training in a two-dimensional field without setting up traditional videoendoscopic camera equipment. It uses patented reflective imaging technology and only requires hand instruments and a suitable tissue medium during use.



Skill Development:

- » Basic laparoscopic skills and dexterity
- » Hand-eye coordination
- » Transferring, dissection and suturing

Features and Benefits:

- » Collapsible for easy storage
- » Portable and affordable
- » No videoendoscopic camera equipment required

Compatible Products:

- » Tissue Model (TSM-10)
- » Complex Tissue Model (TSC-10)
- » Tissue Suture Pad (TSP-10)

Note: Instruments not included