

Torso Trainer

Part #: ST-10

Price: \$625.00 USD

Size: 24" h; 14" w; 18" d

Weight: 8 lbs

This upper abdominal model provides a realistic, durable laparoscopic environment by concealing the procedure from direct view. It accommodates most task and procedure models. The replaceable pneumoperitoneum (ST-20) is shown as an insulated abdomen. Additional replaceable pneumoperitoneum are sold separately.



Skill Development:

- » Laparoscopic training

Features and Benefits:

- » Multi-use trainer
- » Compatible with a variety of Simulab models

Compatible Products:

- » Replaceable pneumoperitoneum (ST-20)
- » Cholecystectomy Model (LCD-10)
- » Key Trainer (SKT-10)