

Thoracentesis Model

Part #: THM-10

Price: \$1,150.00 USD

Size: 12" h; 12" w; 5.5" d

Weight: 6 lbs

This model simulates a partial torso with anatomical landmarks including the scapula, ribs 6 through 10, spinous processes, replaceable skin, pleural cavity, and pleura. The replaceable tissue includes 1/2" of subcutaneous fat. It is designed to be filled with fluid and is self sealing for multiple procedures.

Skill Development:

- » Removal of fluid from pleural cavity
- » Understanding of pleural effusions

Features and Benefits:

- » Anatomically correct with proper landmarks
- » Model allows for full cannulation and the removal of fluid
- » Works in the upright position as well as on its side
- » Model base designed for multiple uses with replaceable tissue
- » Self sealing for multiple procedures

Compatible Products:

- » Replaceable Thoracentesis Tissue and Pleura (THMT-10)

