

TraumaMan® System

Part #: TMPP-30

Price: Request a Quote

Size: 30" h; 22.25" w; 10.5" d (carrying case)

Weight: 22 lbs

The TraumaMan® System, a cutting edge surgical simulator, has been evaluated and approved by the leading Surgical Education Society as an alternative to live non-human models or cadavers for their Surgical Practicum of the leading Trauma Training Course. The TraumaMan® System is an anatomical human body form designed for students to practice several surgical procedures. The System is a simulated human torso with a ventilator and four anatomically correct surgical zones and an ankle base for IV Cutdown. It is designed for replaceable tissue sets that allow each student a "first cut" experience.



Skill Development:

- » Cricothyroidotomy
- » Chest Tube Insertion
- » Pericardiocentesis
- » Needle Decompression
- » Percutaneous Tracheostomy
- » Diagnostic Peritoneal Lavage
- » IV Cutdown

Features and Benefits:

- » Anatomically correct human simulator with all the pertinent landmarks for each procedure
- » The realistic body form skin and replaceable tissue sets have the highest fidelity touch and dissection of any simulated tissue available today
- » TraumaMan is the most flexible training platform available for advanced trauma training
- » Available for purchase, rent or lease

Product Components (purchase):

- » TraumaMan Body Form
- » Replaceable Tissues
- » Ventilator
- » TraumaMan Maintenance Kit
- » Dissected DPL Sample
- » Arterial Fluid (red)
- » Skin Conditioning Powder
- » User's Guide/Instructor's Guide
- » Maintenance DVD
- » Carrying Case