

Complex Tissue Model

Part #: TSC-10

Price: \$75.00 USD

Size: 105mm h; 105mm w; 30mm d

Weight: .5 lbs

The Complex Tissue Model contains skin, subcutaneous fat, fascia, and pre-peritoneal fat. It is an excellent simulation for multi-layer dissecting and suturing.

Skill Development:

- » Dissecting
- » Suturing

Features and Benefits:

- » Realistic
- » Affordable
- » Portable

Compatible Products:

- » Basic Open Surgical Skills (BOSS-10)
- » LapTrainer (LTS-30)

