

# TraumaMan Tribune

## In This Issue

[ACS/APDS Curriculum](#)

[Introducing \[www.traumaman.com\]\(http://www.traumaman.com\)](#)

[Lease-to-Own TraumaMan](#)

[Helpful Hint #5](#)

## Past Issues

[Volume 1 - June 2007](#)

[Volume 2 - Sept 2007](#)

[Volume 3 - Dec 2007](#)

[Volume 4 - March 2008](#)


## Need Expert Assistance?

Contact Simulab's Customer Care Team today!

[Join Our Mailing List!](#)

Welcome to Volume Five of the TraumaMan Tribune. This issue contains a look at TraumaMan and the newly released ACS/APDS Curriculum, an introduction to a new website devoted to TraumaMan, Simulab's new lease-to-own program, and a helpful hint.

## TraumaMan and the ACS/APDS Surgical Skills Curriculum



Simulab is excited about the recently released [ACS/APDS Surgical Skills Curriculum for Residents: Phase I](#) and how we can help meet the [simulation needs for the curriculum](#).

The TraumaMan System can be incorporated into several modules included in the curriculum. These include:

- **Module 1:** Prepping & Draping
- **Module 7:** Surgical Cricothyrotomy and Percutaneous Tracheostomy
- **Module 8:** Chest Tube Insertion

Although TraumaMan does not currently include Percutaneous Tracheostomy, the procedure is one of three upgrades that are currently being incorporated into the TraumaMan System. The upgrade package, which also includes needle decompression and FAST ultrasound, will be available later this year - so stay tuned for more information!

Contact our Customer Care Team at 206-297-1260 or [sales@simulab.com](mailto:sales@simulab.com) for more information.

## Introducing [www.traumaman.com](http://www.traumaman.com)!

**SIMULAB**  
CORPORATION

Realistic, Affordable, and Portable

Home Contact Us

**TraumaMan System**

System Information  
TraumaMan Procedures  
Pricing Options  
Manuals  
Helpful Hints  
TraumaMan Users  
TraumaMan in Action

**TraumaMan Tribune**  
(Click below for past issues)  
Volume 1  
Volume 2  
Volume 3  
Volume 4

Please contact us if you would like to receive the TraumaMan Tribune each quarter.

**TraumaMan® System**  
US Patent #,719,818

For a video on TraumaMan go to our Introduction to the System. (Uses Windows Media Player)

Since **TraumaMan's** debut in 2001, he has both won the approval of the leading surgical societies and become the most widely used surgical simulator in the world. Each year over 20,000 medical professionals are trained using the TraumaMan System and that number keeps growing. TraumaMan was created to meet the needs of advanced trauma surgical skills courses and now has been integrated into a wide variety of curriculum.

The **TraumaMan System** is an anatomical human body form designed for students to practice several advanced surgical procedures, including cricothyroidotomy, chest tube insertion, pericardiocentesis, diagnostic peritoneal lavage, and intravenous cutdown. The **TraumaMan System** consists of simulated human tissue structure made of an elastomeric composition designed specifically for surgical dissection.

To learn more about our **TraumaMan System**, pricing plans, program implementation, or class scheduling, please contact Simulab.

For more information Simulab and other products offered, go to [Simulab's website](http://Simulab's website).

www.traumaman.com

To contact Simulab click here

or Call:  
(206) 297-1260  
(866) 400-1260

To visit Simulab's website go to [www.simulab.com](http://www.simulab.com)

**Trade Exhibits**  
TraumaMan will be at the following shows:

ISEM Annual Meeting in United Kingdom, June 17-21

ACE Clinical Congress in San Francisco, CA, October 12-13

Simulab is pleased to introduce a website devoted entirely to TraumaMan. Visit [www.traumaman.com](http://www.traumaman.com) for information about the system, including procedures, helpful hints, user experiences, acquisition options, and much more.

Please check it out and email [traumaman@simulab.com](mailto:traumaman@simulab.com) if you have ideas or suggestions to improve the site.

### Lease-to-Own TraumaMan

The TraumaMan System is now available on a lease-to-own basis. For customers who are interested in owning the system, but are facing funding or cash flow constraints - you now have another option!

Simulab is offering three different monthly payment plans with varying terms, including 12, 24, 36, 48 or 60 months, to purchase the TraumaMan System.

If you would like more information about this new program, contact Simulab's Customer Care Team at [sales@simulab.com](mailto:sales@simulab.com), 206-297-1260, or 866-400-1260.

## Helpful Hint #5: Chest Tube Insertion



Did you know that each new chest tissue Simulab sends out is accompanied by two new pleura patches? For a realistic experience when performing the chest tube insertion, it is important to discard the used patches and replace them at the same time the chest tissue is changed.

This will maintain the seal around the pleural cavity and allow each student to hear the "pop" of the parietal pleura being penetrated and feel a rush of air from the cavity.

If you would like any additional information about the TraumaMan System or Simulab's other products and services, please call us at 206-297-1260, 866-400-1260, or visit [www.simulab.com](http://www.simulab.com).

Sincerely,



### [Forward email](#)

#### [SafeUnsubscribe®](#)

This email was sent to [sginger@simulab.com](mailto:sginger@simulab.com), by [traumaman@simulab.com](mailto:traumaman@simulab.com)  
[Update Profile/Email Address](#) | Instant removal with [SafeUnsubscribe™](#) | [Privacy Policy](#).

Email Marketing by

