

January 2009 Volume 7

TraumaMan Tribune

In This Issue

[F.A.S.T. Module](#)

[Helpful Hint #7](#)

[Tune up Your System](#)

1st Quarter Special

Visit [Simulab's online storefront](#), enter the code **BOSS9**, and save an additional 10% on the [BOSS Starter Package!](#)

(offer valid until 03.31.09 - online orders only)

Past Issues

[TraumaMan Tribune Vol 1](#)

[TraumaMan Tribune Vol 2](#)

[TraumaMan Tribune Vol 3](#)

[TraumaMan Tribune Vol 4](#)

[TraumaMan Tribune Vol 5](#)

[TraumaMan Tribune Vol 6](#)

[Join Our Mailing List!](#)

Welcome to Volume Seven of the TraumaMan Tribune. This issue contains an overview of the F.A.S.T. Module upgrade for the TraumaMan System and a helpful hint!

Upgrade TraumaMan with the F.A.S.T. Module

The [F.A.S.T. Module for SONORIO](#) is an upgrade for the TraumaMan System that provides an excellent tool to learn F.A.S.T. (Focused Assessment with Sonography for Trauma) exam. This trainer allows the user to learn about:

- Landmarks on the body that are used for identification of window locations for the exam
- Reading and understanding the various normal and abnormal views - F.A.S.T. Module includes five patient image and video libraries for each window (two normal and three abnormal)




With two versions of each image - student and mentor - the user has access to an annotated version of each image. In addition to the mentor version, which provides a standard explanation, the "Instructor" mode allows the instructor to include additional annotations. SONORIO also provides longitudinal views of both the pericardial and pelvic windows and a video for each window of each patient.

For more information, contact Simulab's Customer Care Team: sales@simulab.com or (206) 297-1260.

SONORIO is a trademark of Simulab Corporation.

Helpful Hint #7: Reduce Noise from the Ventilator



There are a couple of steps that can be taken to reduce the noise level of the TraumaMan System ventilator.

The TraumaMan System's ventilator will power two systems at the same time. Even if you are only using one system, always have both hoses attached to the ventilator. Put the hose that is not in use inside the case.

Next, place the ventilator inside the TraumaMan Case and gently close the lid, ensuring that the hoses are not being pinched (folding in the case's latch will prevent this).

Time to Tune up Your System?

Simulab includes a three-year Maintenance Agreement with every TraumaMan System sold. This program includes repair or replacement of the heavily used system components that are vulnerable to excessive incisions, needle sticks, or other damage that can occur during normal training use. If any portion of the system is showing wear and tear, Simulab will refurbish it annually as a part of this program.

In addition, for each 12-month period of the three-year agreement you will receive a new Cricothyroid Cartilage Cover (TMP-1012) free of charge. All repairs and refurbishment under this agreement will be completed at Simulab's facility in Seattle, Washington.

Is it time to get your system refurbished or renew your Maintenance Agreement? Contact our Customer Care Team at 206-297-1260 or sales@simulab.com for instructions.

If you would like any additional information about the TraumaMan System or Simulab's other products and services, please call us at 206-297-1260 or visit [Simulab's website](#).

Sincerely,



[Forward email](#)

 **SafeUnsubscribe®**

This email was sent to sginger@simulab.com by sginger@simulab.com.
[Update Profile/Email Address](#) | Instant removal with [SafeUnsubscribe™](#) | [Privacy Policy](#).

Email Marketing by



Simulab Corporation | 1600 West Armory Way | Seattle | WA | 98119